INSTRUCTIONS FOR USE

TriGUARD 3TM

Cerebral Embolic Protection Device (OUS)

SLB00076 (11)

Contents by language

1	English	1
2	Nederlands	13
3	Deutsch	25
4	Français	37
5	Italiano	49
6	Svenska	61
7	Norsk	73
8	Suomalainen	85
9	Dansk	97
10	Español	109
11	Português	121
12	Čeština	133
13	Slovenščina	145
14	Polish	157
15	Ελληνικά	169
16	Slovak	181

English

INSTRUCTIONS FOR USE

TriGUARD 3™

Cerebral Embolic Protection Device (OUS) SLB00076 (11)


Table of Contents

Table of Contents	1	Physician Training	7
Product Labels	2	Instructions for Use	7
Contents of Package	3	Device Preparation	7
Recommended Supplies	3	Deployment	8
Device Description and Intended Use	4	Aortic Valve and Arch Angiogram	9
Contraindications	5	Trans-Catheter Heart Procedures	9
Warnings	5	Retrieval	10
Adverse Events	6	Disposal	10
Technical Specifications	6	Troubleshooting and Tips	11
How Supplied	7	Warranty	11

Product Labels

The information provided below is an explanation of reference symbols that can be found on product labels.

	Consult Instructions for Use		Contains biological material of animal origin		Non pyrogenic
	Electronic IFU www.keystoneheart.com/eu/ifu		Single sterile barrier system with protective packaging outside.		Keep away from sunlight
	Reference part number		Single sterile barrier system with protective packaging inside		Keep dry
	Lot number		Distributor address		Importer address
	Serial Number		Medical Device		Prescription only device restricted to use by or on the order of a physician
	Sterilized using Ethylene oxide (EO).		Authorized representative in the European community		Temperature limitations
	Caution, consult accompanying documents		Responsible person in the United Kingdom community		Do not re-sterilize
	Do not reuse		Authorized representative in Switzerland		This device has been reviewed and CE Marked
	Manufacture address and date		Do not use if package is damaged		Contains a medicinal substance
	Date of manufacture				Translation
	Use-by-date YYYY-MM-DD				

Contents of Package


Figure 1: Keystone Heart TriGUARD 3™ Cerebral Embolic Protection Device Kit

Recommended Supplies

Not provided in the TriGUARD 3™ Cerebral Embolic Protection Device kit

- Guidewire - Super stiff 0.035", 1cm floppy end, 260 cm long. Ensure curve is added to the straight tip.
- Angio Catheter - 5Fr at least 125cm long, pigtail, or J shaped angio-catheter.
- Pressurized heparinized¹ saline bags

Warning Use only Keystone Heart or procedural recommended components! Use of non-recommended components with the TriGUARD 3 system may result in patient injury, system damage, or loss of efficacy.

¹ Heparinized saline concentration per hospital practice

Device Description and Intended Use

TriGUARD 3 Cerebral Embolic Protection Device Indication

The TriGUARD 3 Cerebral Embolic Protection Device is indicated to minimize the risk of cerebral damage by deflecting embolic debris away from the cerebral circulation during trans-catheter heart procedures.

TriGUARD 3 Cerebral Embolic Protection Device clinical benefit:

The TriGUARD 3 Cerebral Embolic Protection Device is designed to protect the brain from embolic debris.

TriGUARD 3 Cerebral Embolic Protection Device Intended use

Reduce cerebral embolization during transcatheter heart procedures.

TriGUARD 3 Cerebral Embolic Protection Device System Description

The Keystone Heart TriGUARD 3 Cerebral Embolic Protection Device is a cerebral embolic protection device that is designed to provide full coverage of all three major aortic vessels during trans-catheter heart procedures. The TriGUARD 3 Cerebral Embolic Protection Device is designed to minimize the risk of cerebral damage by deflecting embolic debris away from the cerebral circulation during TAVI or other trans-catheter heart procedures.

The TriGUARD 3 Cerebral Embolic Protection Device is introduced through the contralateral femoral artery groin access site and is comprised of an 8Fr delivery system and deflection-filter. It is anatomy independent both for deployment (stability) and for providing coverage since it is designed to accommodate a variety of aortic arch morphologies. The TriGUARD 3 Cerebral Embolic Protection Device is designed to be easily prepared for clinical use, deploy rapidly, self-position, and maintain stability throughout the procedure due to its ability to withstand potential interactions with TAVI delivery systems and accessories.

Target group:

Adult's patients indicated for trans catheter heart procedures.

TriGUARD 3 Cerebral Embolic Protection Device Delivery system

The delivery system (**Figure 2**) is an 8Fr sheath (1) that is compatible with a .035" guidewire. The distal end of the system has a deflection-filter with a nitinol curved hypotube shaft (2) and atraumatic tip (3). The proximal end features a two-part handle (4) with ports for the guidewire (5), pigtail (PT) (6) and heparinized saline flush (7).

The device handle (4) is made of two components, a front part (8) that is connected to the delivery system sheath and rear part (9) that is connected to the hypotube shaft. This design allows relative movements of these two parts to enable the TriGUARD 3 Cerebral Embolic Protection Device to be pulled into the delivery sheath or unsheathed for deployment. The rear part has a luer (5) connected to the hypotube shaft to allow guidewire insertion and flushing (5) of the hypotube shaft and the front part has a place for a tuohy-borst adapter to allow PT insertion (6) and a flushing tube for the 8Fr sheath (7).

TriGUARD 3 Cerebral Embolic Protection Device deflection-filter

The TriGUARD 3 Cerebral Embolic Protection Device deflection-filter (**Figure 3**) is comprised of a frame and mesh. The mesh (**10**) is made of a polymer with a dome shape. The mesh has a nominal pore size of 115 * 145 μm, which allows for adequate blood flow to the cerebral arteries while diverting emboli downstream toward the descending aorta. The frame and mesh are coated with Hydrophilic Heparin coating.

The nitinol frame (11) provides the structural stability of the deflection-filter and is radiopaque for visual confirmation via fluoroscopy. The frame self-expands in the desired location, adapts to the aortic arch anatomy, and provides stability in the aortic arch by inducing radial forces on the aortic arch walls.

The deflection-filter frame ends with a nitinol tail (12). The tail is connected via a nitinol connector to a nitinol curved tube (hypotube Shaft) (2) that has an atraumatic tip at its front end (3). The hypotube shaft crosses underneath the deflection-filter to provide stability and enhanced positioning of the deflection-filter towards the upper part of the aortic arch (passing a stiff guidewire through the hypotube shaft enhances the positioning).

TriGUARD 3 Deflection-Filter

TriGUARD 3 Delivery system

Figure 2: Keystone Heart TriGUARD 3 Cerebral Embolic Protection Device functional parts


Figure 3: Keystone Heart TriGUARD 3 Cerebral Embolic Protection Device deflection-filter

Procedure Overview

The TriGUARD 3 Cerebral Embolic Protection Device is introduced through an 8Fr sheath inserted in one of two femoral arteries at the groin level. Under fluoroscopy, the TriGUARD 3 Cerebral Embolic Protection Device delivery system is inserted over a guide wire and positioned in the ascending arch distal to the innominate artery. Upon deployment from the 8Fr sheath, the deflection-filter unfolds and conforms to the aortic arch anatomy. The frame shape provides apposition to the aortic arch walls to enhance full cerebral vessels coverage (**Figure 4**). Once deployed, the TriGUARD 3 Cerebral Embolic Protection Device can be retrieved into the sheath to allow for repositioning during the procedure or removal at the end of the procedure. The TriGUARD 3 Cerebral Embolic Protection Device is the first system introduced and last system to be removed after the index procedure is completed.

The TriGUARD 3 is supplied sterile and should not be re-sterilized.

Contraindications

DO NOT USE IN THE FOLLOWING CIRCUMSTANCES:

- Hypercoagulable states that cannot be corrected by additional periprocedural heparin.
- Renal failure with plasma creatinine > 4 mg/dl.
- Hepatic failure.
- Patients with allergy to nitinol or heparin
- Patients with history of heparin-induced thrombocytopenia


Figure 4: TriGUARD 3 Cerebral Embolic Protection Device position in the aortic arch

Warnings

- Only physicians trained in the use of the device should use it.
- The appropriate antiplatelet and anticoagulation therapy should be administered pre- and post-procedure to minimize the risk of embolism and thrombosis.
- Porcine material as porcine allergens other than porcine heparin could be present
- Failure to follow recommended device preparation and use of dry pad to wipe the filter may damage the hydrophilic coatings and potentially cause serious injuries to patients.
- Introduce and advance devices slowly to prevent air embolism or trauma to the vasculature.
- Do not attempt to move the handle without observing the resultant deflection-filter response.
- Visualization of device position is essential for proper deployment. In procedures where visualization of the device via fluoroscopy is impeded by patient obesity or other reasons, do not attempt to deploy the device in the aortic arch.
- Ensure visualization of the device when attempting to cross the aortic arch with any device. Uncontrolled interaction between the TriGUARD 3 Cerebral Embolic Protection Device and any other device may lead to undesirable movement of the deflection-filter, resulting in incomplete coverage of the ascending arteries.
- To avoid damage to the device or injury to the patient, do not pull or push excessively on the device handle or the Introducer Sheath.
- Rotation of the device handle may result in rotation or flipping of the device. Do not attempt to rotate or torque the device using any accessory tool or the device handle.
- TriGUARD 3 is a single use device. Do not reuse or re-sterilize. Reusing the device may impose risk of inter-patient contamination, improper cleaning, and compromised

performance of the device.

- In case of redeployment, first pull back the device into the sheath and then advance the sheath over the wire to the desired location. Once in the desired location, deploy per the instructions for use.
- The safety and effectiveness of the device has not been evaluated in patient populations with the following conditions:
 - Variant angina pectoris, unstable angina or recent acute coronary syndrome including myocardial infarction (in the past three months).
 - Stroke, TIA - Transient Ischemic Attack (in the past 6 months).

Adverse Events

There are risks associated with any endovascular intervention and complications may develop. The following anticipated events have been identified as possible complications of trans-catheter procedures and these and others may be associated with the device:

- Acute cardiovascular surgery (need for)
- Acute coronary artery occlusion
- Acute myocardial infarction
- Acute neurological events such as: Stroke, transient ischemic attack (TIA), encephalopathy
- Allergic reaction to contrast, antiplatelet therapy, or device component materials
- Bowel ischemia
- Embolism (air, tissue, device, or thrombus)
- Fever
- Femoral nerve damage
- Fluoroscopy related harm
- Hemodynamic changes
- Angina pectoris
- Anesthesia reactions
- Aneurysm or pseudoaneurysm
- Arteriovenous fistula
- Ascending or descending aorta trauma
- Atrial or ventricular arrhythmias or fibrillation, Heart Palpitations (sustained requiring therapy)
- Bleeding complications such as hematoma and hemorrhage
- Bleeding at access site
- Blood loss requiring transfusion
- Coronary artery or other vascular injury, dissection, or perforation which may need repair
- Blue toe syndrome or blue discoloration of a toe
- Hypertension or hypotension (sustained requiring therapy)
- Infection, including endocarditis and septicemia
- Pain (at femoral puncture site, abdominal, back, or other)
- Percutaneous coronary intervention (need for)
- Peripheral ischemia, peripheral nerve damage

- Hypotension (systolic blood pressure of below 90mm Hg).
- Active peptic disease or history of upper GI bleeding.
- Spastic bronchitis, Chronic Obstructive Lung Disease, Asthma.
- Complex ventricular arrhythmia or history thereof.
- Major psychiatric disorder in the present or past.
- Bleeding diathesis such as hemophilia, ITP, aplastic anemia, TTP etc.
- Diabetes mellitus with peripheral vascular and/or neurologic changes.
- Any proliferative disease with patient life-span less than 6 months.
- Pregnancy.
- Pediatric use.
- Pulmonary edema
- Pyrogenic reaction
- Renal complications, injury, or failure
- Unstable angina
- Vascular complications which may require vessel repair
- Vessel spasm (sustained, not responding to therapy)

In addition to the risks listed above, the potential risks specifically associated with the Keystone Heart TriGUARD 3 Cerebral Embolic Protection Device procedure includes, but may not be limited to, the following:

- Dislodgement or migration of the TriGUARD 3 Cerebral Embolic Protection Device or its delivery system, due to passage of other instrumentation, e.g.: balloon, stent, catheter, wire.
- Femoral bleeding at the access site.
- Local trauma to the aortic wall due to device movement.
- Any serious incident occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State.

Technical Specifications

Deflection-filter	Frame width	74 mm
	Frame length	98 mm
Delivery hypotube shaft	Inner diameter	Accommodates 0.035" Guide-wire
	Total length	127.5 cm
8Fr Sheath	Inner diameter	8 Fr
	Effective length	76 cm
Device effective length	Strain relief to atraumatic tip (during over the wire advancement)	78cm

SHELF LIFE: 12 months

How Supplied

Packaging

The kit is pre-assembled and packed in an inner blister tray with a blister cover. The blister is placed in a sealed Tyvek pouch. The set is placed in a cardboard box. The entire set is Ethylene Oxide (EtO) sterilized. Only the blister Tray should be placed in the sterile field. Inspect packaging prior to use. Do not use if there is damage to the packaging, the package is open, or if there are any other defects. Ensure appropriate labeling and that the device is used prior to device expiration.

Storage

Store at room temperature in a dry, dark (away from sunlight) place.

Physician Training

Only physicians trained in percutaneous intravascular techniques and procedures should use the Keystone Heart TriGUARD 3 Cerebral Embolic Protection Device System.

Keystone Heart TriGUARD 3 Cerebral Embolic Protection Device System training is required and may include on-site training and proctoring of initial cases.

Instructions for Use

Device Preparation

1. Inspect the package sealing and verify product sterility and integrity. An opened or damaged item should not be used and should be returned to Keystone Heart.
2. Open the TriGUARD 3 Cerebral Embolic Protection Device carton box.
3. Open the TriGUARD 3 Cerebral Embolic Protection Device sealed Tyvek pouch and place the sterile blister tray in the sterile zone.
4. Remove the blister tray cover.
5. Tighten the handle nut to the front part of the handle to a full closure (**Figure 5 - 5A**).

6. Fill in the flushing basin with saline (or heparinized saline) until the TriGUARD 3 deflection-filter is fully immersed.
7. Flush the hypotube shaft with saline (or heparinized saline) through the luer located at the rear part of the device handle (**Figure 5 - 5B**).
8. Flush the delivery system sheath with saline (or heparinized saline) through the flushing tube located at the front part of the device handle (**Figure 5 - 5C**).
9. Once saline (or heparinized saline) is dripping out of the Tuohy Borst, screw it to full closure (**Figure 5 - 5D**).
10. Immerse the TriGUARD 3 deflection-filter in the flushing basin for approximately 1-minute in saline (or heparinized saline) solution to hydrate the heparinized hydrophilic coating.
11. After 1 minute of immersion, no air bubbles should be visible. Gently tap on the deflection filter to remove remaining air bubbles.
12. While maintaining the handle orientation, pullback the rear part of the handle while holding the **front part** stationary until the Deflection-Filter is fully crimped into the delivery system.
13. Ensure that the atraumatic tip is not fully inside the sheath or extending beyond the tip of the sheath. If a gap is noted between the atraumatic tip and the sheath, address by advancing the sheath over the atraumatic tip.
14. Flush the 8Fr sheath via the flushing tube with saline (or heparinized saline), while TriGUARD 3 Cerebral Embolic Protection Device tip is immersed completely, until no bubbles are released.
15. Pull back the device delivery system until delivery sheath is totally out of the protecting sleeve.

Figure 5: Flushing the hypotube shaft via the Luer and the 8Fr sheath via the flushing tube

Note Handle orientation should be kept while crimping the deflection filter into the delivery system by: Tuohy and flushing tube are facing toward the left side, as indicated in **Figure 5**. The delivery system tip is inserted into the crimper.

Warning Avoid wiping the TriGUARD 3 Cerebral Embolic Protection with dry gauze as this may damage the device coating.

Deployment

Warning Perform all steps under fluoroscopic guidance. Do not make any movements of the delivery system or components without adequate visualization. All steps that reference delivery system movements assume fluoroscopic guidance!

Warning Failure to maintain an ACT level >250 may increase the risk of thrombus formation on the device and in the 8Fr delivery sheath.

1. Maintain an ACT>250 seconds throughout device deployment and dwell time, in line with routine hospital practice.
2. Insert a guidewire in the ascending Aorta in the vicinity of native aortic annulus.
3. In a wide LAO projection, perform an angiogram to demonstrate the anatomy of the aortic arch and cerebral vessels. If possible, superimpose the aortic arch image on top of the live fluoroscopic image for reference purposes. It is recommended to obtain the best anatomical view of the arch and cerebral vessels. Remove angiogram catheter once completed.

Note Before insertion of the delivery system into the introducer/groin, it is recommended to inject Heparinized saline to the delivery system through the flushing tube until the tip of the delivery system drips

4. Before insertion of the delivery system into the introducer groin, inject HS through the flushing tube until the tip of the delivery sheath drips (use 5-10cc syringe).
5. Advance the TriGUARD 3 Cerebral Embolic Protection Device delivery system (loaded with the TriGUARD 3 Cerebral Embolic Protection Device), over the guidewire, to approximately 4 cm distal to the Innominate branch in the ascending aorta (based on the initial angiogram) using fluoroscopic guidance in the same reference plane as the angiogram (**Figure 6**).


Figure 6: The delivery sheath is positioned 4 cm distal to the Innominate branch in the ascending aorta

Note Insertion of the TriGUARD 3 Cerebral Embolic Protection Device into the femoral access should be made by advancing the 8Fr sheath only.

Warning Avoid manually rotating the TriGUARD 3 Cerebral Embolic Protection Device deflection-filter either while inside the delivery sheath or at the deployment site.

Warning Advance the TriGUARD 3 Cerebral Embolic Protection Device, over the wire, while avoiding torquing.

Warning The handle to rotate freely according to vessel trajectory during the insertion of the delivery system

6. Pull the guidewire back to the descending aorta to confirm the hypotube shaft is under the mesh when deployed and to allow the mesh to fully expand across the aortic arch
7. Hold the rear part of the handle stationary and slowly pull back the front part of the handle until the TriGUARD 3 Cerebral Embolic Protection Device deflection-filter is fully deployed from the delivery sheath (**Figure 7**). Re advance the guidewire to the vicinity of the native aortic annulus.
8. Connect a continuous, pressurized heparinized saline bag to the delivery-sheath flushing tube located in the TriGUARD 3 Cerebral Embolic Protection Device handle. Another option is to perform a heparinized saline flush (every 20 min) through the flushing tube to maintain the sheath free from blood clots.
9. Maintain forward pressure on the guide-wire of the TriGUARD 3 Cerebral Embolic Protection Device at the annulus. This will provide support for the mesh and keep the atraumatic tip on the outer curvature of the arch.


Figure 7: The TriGUARD 3 Cerebral Embolic Protection Device is deployed in the Aortic arch

Note Following deployment, pull and push the GW in order to verify that the hypotube is underneath the deflection filter and to allow the deflection filter to self-position.

Note User may control the hypotube shaft curvature by pulling/pushing the GW in the hypotube: Pushing the guidewire will push the hypotube shaft toward the upper curvature of the arch providing better support for the deflection filter and moving the tip away from mid-stream.

Warning Leave sufficient guide-wire distal to the atraumatic tip to avoid any potential damage to the vessel wall from the guidewire or tip.

Note If the guidewire is pulled out from the TriGUARD 3 Cerebral Embolic Protection Device delivery system, it is recommended to connect the Luer, at the rear part of device's handle, to a continuous pressurized dripping system, to maintain the hypotube shaft free from blood clots.

Aortic Valve and Arch Angiogram

Note It is recommended to perform an angiogram using a pig-tail angiography catheter.

1. Advance the pigtail catheter through the rotating tuohy-borst adapter at the front part of the TriGUARD 3 Cerebral Embolic Protection Device handle. Use fluoroscopy to assure its position underneath the deflection-filter by verifying that it does not alter the position of the deflection-filter during tracking.

Note It is best to introduce the pigtail with the curled end leading under the deflection-filter Mesh. It is recommended to use at least two views to ensure crossing underneath the deflection-filter.


Figure 8: Deflection-filter optimal coverage and apposition

2. After passing the deflection-filter, advance the pigtail catheter into ascending aorta. It is recommended to perform an angiogram to evaluate device position in the aortic arch. If needed, fine adjustments of device position may be made for optimal coverage.

Trans-Catheter Heart Procedures

The TriGUARD 3 may be used in accordance with the indications for use in transcatheter heart procedures including but not limited to:

- **Trans-catheter Aortic Valve Implantation (TAVI)**
- **Left Atrial Appendage Closure (LAAC)**
- **Atrial Fibrillation Ablation**
- **Transcatheter Mitral Valve Repair/Replacement (TMVR)**

Follow all index procedural related steps per proper device instructions for use.

Note For best results, the TriGUARD 3 Cerebral Embolic Protection Device should be First-in Last Out (FILO) to be sure that full vessel coverage and protection is maintained throughout the index procedure.

When an attempt is made to pass the aortic arch with any additional devices or accessories (e.g. guide-wire, TAVI, balloon) make sure that:

- The deflection-filter covers the aortic branches with its front and backend opposing the aortic walls. (**Figure 8**).
 - The hypotube shaft with the atraumatic tip at its distal part is under and opposing the center of the deflection-filter.
1. Advance the guidewire (ipsilateral access) using fluoroscopy to verify that the guidewire is below the deflection-filter and does not alter its position.
 2. At physician discretion, complete the main index procedure.

Note When advancing instrumentation (balloons, guidewires, valve delivery systems or other devices) around the arch and under the deflection-filter, counterpull the guidewire while pushing the delivery system forward. This will prevent tension build-up and will allow smoother passage of the instrumentation under the TriGUARD 3.

Warning Always visualize TriGUARD 3 during the entire procedure to ensure the TAVR catheter does not get tangled with deflection filter potentially causing dislodgement. **Uncontrolled interaction between the TriGUARD 3 Cerebral Embolic Protection Device and any other device may lead to inadvertent movement of the deflection-filter.**

Warning Trying to manipulate/ crossing underneath the TriGUARD 3 deflection filter with a bare wire is forbidden.

Warning When TAVI or any other delivery system is retrieved (pulled back), make sure it does not entangle with the TriGUARD 3 Cerebral Embolic Protection Device deflection-filter and pulls it back.

Retrieval

1. Remove the trans-catheter devices used during the procedure and the pigtail catheter from the TriGUARD 3 Cerebral Embolic Protection Device delivery system.

Note User shall retrieve the pigtail catheter over a guide wire to avoid entanglement with the deflection-filter.

2. Pull back the handle from its rear part, while holding the front handle front part stationary, until the TriGUARD 3 deflection filter is fully crimped into the delivery sheath (**Figure 9**).
3. Pull back the delivery sheath, with the crimped TriGUARD 3 Until the delivery sheath is fully removed from the patient's body.
4. Close femoral access point.

Disposal

The TriGUARD 3 Device must be disposed of in accordance with local biohazard waste disposal and hospital procedures.


Figure 9: The TriGUARD 3 Cerebral Embolic Protection Device retrieval

Troubleshooting and Tips

Note Events described below may take place during the procedure and can be handled using the following trouble shooting tips.

1. TriGUARD 3 Cerebral Embolic Protection Device malpositioned during deployment or migration in the transverse arch.

Retrieve the deflection-filter into the delivery sheath by pulling back the rear part of the TriGUARD 3 Cerebral Embolic Protection Device deflection-filter handle while holding its front part stationary until the TriGUARD 3 Cerebral Embolic Protection Device deflection-filter is fully crimped into the delivery sheath and the device tip meets the sheath tip. Reposition the delivery sheath, over the guide wire, to the correct location and then redeploy as described above.

Note Only one additional attempt of complete redeployment of the same device is permitted.

Note Make sure both the introducer sheath and the delivery sheath are secured in position.

2. The position of the guidewire or pigtail catheter underneath the TriGUARD 3 Cerebral Embolic Protection Device is not clear.

Using fluoroscopy (RAO view), slightly pull and push guidewire or pigtail catheter while observing its peak location. If this manipulation causes the guidewire or pigtail catheter to be in the same plane as the hypotube shaft, then guidewire or pigtail catheter is underneath the TriGUARD 3. If it is NOT in the same plane, remove the guidewire or pigtail and re-navigate it below the deflection-filter Mesh.

3. Deflection-filter is out of position during maneuvering of other devices in the aortic arch.

Check if the introducer sheath and the delivery sheath are secured in position and slowly advance index procedure hardware.

If deflection-filter fails to return to its position, remove it before continuing with the index procedure or

At physician discretion, it may be attempted to retract the index procedure hardware to the descending aorta and reposition the deflection-filter.

4. Interaction between the TriGUARD 3 and TAVI delivery system.

Slightly pull back (~1 cm) the TAVI delivery system and the TAVI wire which will disengage the TAVI tip from the TriGUARD 3 device frame. Then push the TAVI delivery system forward while holding stationary the TAVI wire which will lower the tip of the TAVI delivery system to allow free passage underneath the TriGUARD 3 device frame.

Warranty

Keystone Heart warrants that reasonable care has been used in the design and manufacture of this instrument. The foregoing warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied, by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.

Handling, storage, cleaning and sterilization of this instrument, as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond Keystone Heart's control directly affect the instrument, its performance and the results obtained from its use.

Keystone Heart's obligation under this limited warranty is strictly limited to the replacement of the instrument. In no event will Keystone Heart be liable for any indirect, incidental or consequential loss, damage or expense directly or indirectly arising from the use of the instrument. Keystone Heart neither assumes nor authorizes any other person to assume for it any other or additional liability or responsibility in connection with the instrument. Keystone Heart assumes no liability with respect to any re-use, re-processing or re-sterilization of the instrument, and disclaims all warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.

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Switzerland Authorized Representative

Obelis Swiss GmbH
Registered Address
Ruessenstrass 12, 6340
Baar/ZG Switzerland
Tel: +(04)15441.526
Fax: +(04)15441.527
E-mail: info@obelis.ch

European Authorized Representative

Obelis s.a
Bd. Général Wahis 53
1030 Brussels, BELGIUM
Tel: +(32)27325.954
Fax: +(32)27326.003
E-mail: mail@obelis.net

Responsible Person

Obelis UK LTD
Sanford Gate
East Point Business Park
Oxford OX4 6LB
United Kingdom
Tel: +(44)1865910.145
E-mail: info@Obelis.co.uk

Swiss Importer
MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

EU Importer
MedEnvoy
Prinses Margrietplantsoen 33
Suite 123
2595 AM The Hague
The Netherlands

UK Importer
MedEnvoy UK Limited
85,Great Portland Street
First Floor
London, W1W7LT
United Kingdom

Keystone Heart LTD
Halamish St. 15
3088900 Caesarea
Israel
Tel: +(972)4615.8000
Fax: +(972)4615.8099
www.keystoneheart.com

Translation Company
Synergy Language Services Ltd.
Parkhill Business Centre Annexe
Walton Road
Wetherby West Yorkshire
LS22 5DZ
United Kingdom
Tel: +(44)01937.547.069
E-mail: sales@synergylanguageservices.com